



Product Specification
Screw cap tube, 10 ml, (LxØ): 92 x 15.3 mm, PP,
with print

Page 1



Product description

Order number	62.610.201
Product description	Screw cap tube, working volume: 10 ml, (LxØ): 92 x 15.3 mm, material: PP, round base, transparent, screw cap, transparent, cap assembled, with print, label/print: white, with graduations, sterile, 100 piece(s)/bag

Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	round base
Volume of work	10 ml
Centrifugation max (RCF)	3000 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Screw cap tube, 10 ml, (LxØ): 92 x 15.3 mm, PP,
with print**

Page 2

Size

Diameter	15.3 mm
Length of product	92 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	transparent

Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0.0306 cbm
Weight of product	0.0052 kg
Weight of case	3.13 kg
EAN of inner box	4038917114783
EAN case	4038917072137

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

