



Product description

Order number	63.504.027
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, PP, transparent, cap assembled, with plastic label, label/print: transparent/green, "CSF/Liquor", with graduations, cap assembled, natural, sterile, 1 piece(s)/blister

Product characteristics

Label/ Print	with plastic label
Colour of print/label	transparent/green
Graduation	yes
Closure type	screw cap
Example of use	CSF, collection, transport and storage of CSF/liquor
Base shape	round base
Volume of work	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Diameter	13 mm
Length of product	75 mm
Length including cap	83 mm
Length excluding cap	75 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	13500
Depth of box	351 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.0039 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Screw cap tube, 5 ml, round base, transparent

Weight of case	3.56 kg
EAN of inner box	4038917408561
EAN case	4038917408554

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

