



## Product Specification

**Tube with false bottom CSF, 2.5 ml, (LxØ): 75 x 13 mm, round base, PP, Low Binding, cap enclosed, sterile, 1 piece(s)/blister**

**Page 1**



### Product description

Order number	63.614.625
Product description	Tube with false bottom CSF, nominal volume: 2.5 ml, (LxØ): 75 x 13 mm, round base, transparent, material: PP, Low Binding, with plastic label, label/print: transparent/blue, "CSF", fill mark, cap enclosed, blue, sterile, 1 piece(s)/blister, 100 piece(s)/carton

### Product characteristics

Label/ Print	with plastic label
Colour of print/label	transparent/blue
Graduation	fill mark
Light protection	no
Application	collection and storage of diagnostic specimens
Example of use	Alzheimer's disease (AD) diagnostics
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	20

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Tube with false bottom CSF, 2.5 ml, (LxØ): 75 x 13 mm, round base, PP, Low Binding, cap enclosed, sterile, 1 piece(s)/blister**

**Page 2**

Centrifugation time 10 Min  
Shape of inner tube conical

### Size

Sample volume 2.5 ml  
Diameter 13 mm  
Length of product 75 mm

### Material & colours

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap blue  
Tube material Polypropylene (PP)

### Purity & certification

Product category In-vitro diagnostic (IVD)  
CE certified CE - manufacturer's self-declaration  
Purity standard sterile  
Sterilisation Gas sterilisation  
Batched yes

### Packaging

Minimum order qty. 100  
Type of smallest subpackaging blister  
Piece(s) / inner box 10  
Piece (s) / outer case 100  
Piece(s) / pallet 12100  
Depth of case 338 mm  
Width of case 230 mm  
Height of case 153 mm  
Case volume 0.0119 cbm  
Weight of product 0.0038 kg  
Weight of case 0.82 kg  
EAN of inner box 4038917896641

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Tube with false bottom CSF, 2.5 ml, (LxØ): 75 x 13 mm, round base, PP, Low Binding, cap enclosed, sterile, 1 piece(s)/blister**

**Page 3**

EAN case

4038917896658

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**