



Product description

| | |
|---------------------|---|
| Order number | 95.1000.100 |
| Product description | Tube segment opener, for bed side test cards, material: HD-PE, white, 200 piece(s)/case |

Product characteristics

| | |
|-------------|-------------------------|
| Design | for bed side test cards |
| Application | Blood bank |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | High Density Polyethylene (HD-PE) |
| Colour of product | white |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Tube segment opener, for bed side test cards,
HD-PE, white

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 200 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 242000 |
| Depth of box | 139 mm |
| Width of box | 108 mm |
| Height of box | 66 mm |
| Depth of case | 338 mm |
| Width of case | 228 mm |
| Height of case | 160 mm |
| Case volume | 0.0123 cbm |
| Weight of product | 0.0018 kg |
| Weight of case | 4.48 kg |
| EAN of inner box | 4038917124607 |
| EAN case | 4038917028585 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

