



Product Specification
Tube segment opener, for tubes with inner Ø 8.2 mm, HD-PE, white

Page 1



Product description

Order number	95.1000.800
Product description	Tube segment opener, for tubes with inner Ø 8.2 mm, material: HD-PE, white, 500 piece(s)/case

Product characteristics

Application	Blood bank
Suitable for	tubes with inner Ø 8.2 mm

Material & colours

Product material	High Density Polyethylene (HD-PE)
Colour of product	white

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Tube segment opener, for tubes with inner Ø
8.2 mm, HD-PE, white

Page 2

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	216000
Depth of box	191 mm
Width of box	129 mm
Height of box	129 mm
Depth of case	398 mm
Width of case	268 mm
Height of case	150 mm
Case volume	0.016 cbm
Weight of product	0.0015 kg
Weight of case	3.6 kg
EAN of inner box	4038917124645
EAN case	4038917053273

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT