



Product description

Order number	93.1650.010
Product description	Insertion plate, Version 1, material: PP, 1 piece(s)/case

Product characteristics

Suitable for	Safety Tray (93.1650) & e.g. 1x Multi-Safe twin plus 500 (77.3884.005/.105), 1x rack with handle (93.851 - 93.851.1xx), 2x dispenser boxes (95.813.001), 3x multi-purpose cup (75.1355), 1x disinfection spray
--------------	--

Size

Width of product	272 mm
Length of product	420 mm

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1440
Case volume	0.0009 cbm
Weight of product	0.152 kg
Weight of case	0.31 kg
EAN of inner box	4038917132695
EAN case	4038917076715

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

