



Product description

Order number 93.1650.030
Product description Insertion plate, Version 3, material: PP, 1 piece(s)/case

Product characteristics

Suitable for Safety Tray (93.1650) & e.g. 1x Multi-Safe twin plus 500 (77.3884.005/.105), 2x dispenser boxes (95.813.001), 8x multi-purpose cup (75.1355), 1x disinfection spray

Size

Width of product 272 mm
Length of product 420 mm

Material & colours

Product material Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 1 |
| Piece(s) / pallet | 1440 |
| Case volume | 0.0009 cbm |
| Weight of product | 0.172 kg |
| Weight of case | 0.33 kg |
| EAN of inner box | 4038917132718 |
| EAN case | 4038917076739 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

