



Product description

Order number	95.813.001
Product description	Dispenser box, material: PS, lid assembled

Product characteristics

Base shape	flat base
------------	-----------

Size

Diameter of base	81 mm
Diameter of lid	83 mm
Diameter of opening	81 mm
Length of product	110 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	white
Closure material	Polystyrene (PS)
Colour of cap	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 35
Piece(s) / inner box 35
Piece (s) / outer case 35
Piece(s) / pallet 1260
Depth of case 398 mm
Width of case 398 mm
Height of case 300 mm
Case volume 0.0475 cbm
Weight of product 0.0623 kg
Weight of case 2.77 kg
EAN of inner box 4038917121644
EAN case 4038917074797

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

