




## Product Specification

Page 1

Round tube holder for round tubes and combinable with wall holder euroMatic® (95.963.001), wall holder twin plus (95.963.007) and wall holder vario 2000 (95.963.020)



### Product description

Order number	95.963.006
Product description	Round tube holder, for round tubes and combinable with wall holder euroMatic® (95.963.001), wall holder twin plus (95.963.007) and wall holder vario 2000 (95.963.020), material: POM, (LxWxH): 130 x 40 x 88 mm

### Product characteristics

Suitable for	round tubes and combinable with wall holder euroMatic® (95.963.001), wall holder twin plus (95.963.007) and wall holder vario 2000 (95.963.020)
--------------	---

### Size

Width of product	40 mm
Height of product	88 mm
Length of product	130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

Page 2

Round tube holder for round tubes and combinable with wall holder euroMatic® (95.963.001), wall holder twin plus (95.963.007) and wall holder vario 2000 (95.963.020)

### Material & colours

Product material	Polyoxymethylene (POM)
Colour of product	silver/white

### Purity & certification

Product category	no medical device   IVD
Batched	no

### Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	1680
Depth of case	288 mm
Width of case	218 mm
Height of case	150 mm
Case volume	0.0094 cbm
Weight of product	0.4003 kg
Weight of case	4.23 kg
EAN of inner box	4038917364300
EAN case	4038917364294

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**