



### **Product description**

Order number	95.963.007
Product description	Wall holder, for Multi-Safe twin plus, material: stainless steel

### **Product characteristics**

Example of use	disposal of pointed and sharp objects
Suitable for	Multi-Safe twin plus

### **Size**

Width of product	186 mm
Length of product	165 mm

### **Material & colours**

Product material	stainless steel
------------------	-----------------

### **Purity & certification**

Product category	no medical device   IVD
Latex-free	yes
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	20
Piece(s) / pallet	2880
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.246 kg
Weight of case	5.24 kg
EAN of inner box	4038917900997
EAN case	4038917901000

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

