




Product description

Order number	95.963
Product description	Wall holder, for Multi-Safe twin plus, euroMatic® , with adhesive/flat adapter assembled, material: stainless steel

Product characteristics

Design	adhesive/flat adapter assembled
Suitable for	Multi-Safe twin plus, euroMatic®

Size

Width of product	51 mm
Length of product	139 mm

Material & colours

Product material	stainless steel
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 30
Piece(s) / pallet 2880
Depth of case 348 mm
Width of case 278 mm
Height of case 140 mm
Case volume 0.0135 cbm
Weight of product 0.1801 kg
Weight of case 5.81 kg
EAN of inner box 4038917084697
EAN case 4038917003469

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

