



**Product description**

Order number	77.3884.120
Product description	Disposal container, Multi-Safe twin plus, 2,000 ml, yellow, with plastic label, label/print: yellow/red, biohazard labelling, material: PP, application e.g.: disposal of pointed and sharp objects, opening: oval, height: 161 mm, lid (LxW): 179 x 176 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 80 piece(s)/carton

**Product characteristics**

Label/ Print	with plastic label
Colour of print/label	yellow/red
Example of use	disposal of pointed and sharp objects
Shape of opening	oval
Removal aid - turning aid	Removal and rotation aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Max. volume	2000 ml
Diameter of base	112 mm
Height of product	161 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Closure material	Polypropylene (PP)
Colour of cap	red
Colour of lid	red

### **Purity & certification**

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	not regulated
CE certified	no CE mark
Batched	no

### **Packaging**

Minimum order qty.	80
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	80
Piece(s) / pallet	960
Depth of case	800 mm
Width of case	400 mm
Height of case	380 mm
Case volume	0.1216 cbm
Weight of product	0.194 kg
Weight of case	17.17 kg
EAN of inner box	4038917903448
EAN case	4038917903455

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

