



### **Product description**

Order number	77.3890.003
Product description	Disposal container, Multi-Safe mini, 300 ml, yellow, material: PP, opening: round, height: 160 mm, flat cap, lid (LxW): 81 x 45 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 135 piece(s)/carton
UN identification	yes

### **Product characteristics**

Cap	flat cap
Shape of opening	round
Removal aid - turning aid	Removal aid

### **Size**

Max. volume	300 ml
Diameter of base	80 mm x 45 mm
Diameter of opening	27 mm
Height of product	160 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	red

### **Purity & certification**

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	135
Type of smallest subpackaging	case
Piece(s) / inner box	135
Piece (s) / outer case	135
Piece(s) / pallet	2700
Depth of case	400 mm
Width of case	310 mm
Height of case	600 mm
Case volume	0.0744 cbm
Weight of product	0.0694 kg
Weight of case	10 kg
EAN of inner box	4038917330947
EAN case	4038917330930

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

