



Product description

Order number	77.3897.030
Product description	Disposal container, Multi-Safe euroMatic® , 3,000 ml, yellow, square, with plastic label, label/print: yellow/red, material: PP, opening: oval, height: 185 mm, Automatic lock, lid (LxW): 164 x 164 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 55 piece(s)/carton
UN identification	yes

Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Cap	Automatic lock
Shape of opening	oval
Removal aid - turning aid	Removal aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Max. volume	3000 ml
Diameter of base	132 mm x 132 mm
Diameter of opening	110 mm x 60 mm
Height of product	185 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	red

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	55
Type of smallest subpackaging	case
Piece(s) / inner box	55
Piece (s) / outer case	55
Piece(s) / pallet	495
Depth of case	392 mm
Width of case	514 mm
Height of case	802 mm
Case volume	0.1616 cbm
Weight of product	0.2835 kg
Weight of case	17 kg
EAN of inner box	4038917332897
EAN case	4038917332880

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

