



Product description

Order number	77.3897.270
Product description	Disposal container, Multi-Safe euroMatic® , 7,000 ml, yellow, with plastic label, label/print: yellow/red, biohazard labelling, material: PP, opening: oval, height: 390 mm, flat cap, lid (LxW): 164 x 164 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 34 piece(s)/carton
UN identification	yes

Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Cap	flat cap
Shape of opening	oval
Removal aid - turning aid	Removal aid

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

Max. volume	7000 ml
Diameter of base	132 mm x 132 mm
Diameter of opening	110 mm x 60 mm
Height of product	390 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	red

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Minimum order qty.	34
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	34
Piece(s) / pallet	306
Depth of box	392 mm
Width of box	514 mm
Height of box	802 mm
Depth of case	802 mm
Width of case	392 mm
Height of case	514 mm
Case volume	0.1616 cbm
Weight of product	0.5228 kg
Weight of case	19 kg
EAN of inner box	4038917963046
EAN case	4038917963053

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

