



Product Specification
Disposal container, Multi-Safe euroMatic® , 6,000 ml, biohazard labelling

Page 1

Product description

Order number	77.3897.660
Product description	Disposal container, Multi-Safe euroMatic® , 6,000 ml, grey, square, with plastic label, label/print: yellow/red, biohazard labelling, material: recycled PP, opening: oval, height: 330 mm, Automatic lock, lid (LxW): 164 x 164 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 40 piece(s)/carton
UN identification	yes

Product characteristics

Label/ Print	with plastic label
Colour of print/label	yellow/red
Cap	Automatic lock
Shape of opening	oval
Removal aid - turning aid	Removal aid

Size

Max. volume	6000 ml
Diameter of base	132 mm x 132 mm
Diameter of opening	110 mm x 60 mm
Height of product	330 mm

Material & colours

Product material	recycled Polypropylene (recycled PP)
Colour of product	grey
Colour of cap	red
Colour of lid	red

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	not regulated
CE certified	no CE mark
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Disposal container, Multi-Safe euroMatic® , 6,000 ml, biohazard labelling

Page 2

Packaging

Minimum order qty.	40
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	40
Piece(s) / pallet	360
Depth of case	802 mm
Width of case	392 mm
Height of case	514 mm
Case volume	0.1616 cbm
Weight of product	0.4398 kg
Weight of case	19 kg
EAN of inner box	4066916004965
EAN case	4066916004972

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-03-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT