



Product description

| | |
|---------------------|---|
| Order number | 77.3898.111 |
| Product description | Disposal container, Multi-Safe medi, 11 l, yellow, round, material: PP, with assembled bracket, opening: round, height: 220 mm, flat cap, Ø of cap: 290 mm, cap: red, TRBA 250, DIN EN ISO 23907, 250 piece(s)/carton |
| UN identification | no |

Product characteristics

| | |
|---------------------------|------------------------|
| Cap | flat cap |
| Design | with assembled bracket |
| Shape of opening | round |
| Removal aid - turning aid | Removal aid |

Size

| | |
|---------------------|--------|
| Max. volume | 11 l |
| Diameter of base | 257 mm |
| Diameter of lid | 290 mm |
| Diameter of opening | 88 mm |
| Height of product | 220 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | yellow |
| Colour of cap | red |
| Colour of lid | yellow |

Purity & certification

| | |
|---------------------------|----------------------------|
| Satisfies the requirement | TRBA 250, DIN EN ISO 23907 |
| Product category | not regulated |
| CE certified | no CE mark |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 250 |
| Type of smallest subpackaging | pallet |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 250 |
| Piece(s) / pallet | 250 |
| Depth of case | 1192 mm |
| Width of case | 792 mm |
| Height of case | 1090 mm |
| Case volume | 1.029 cbm |
| Weight of product | 0.6135 kg |
| Weight of case | 163.32 kg |
| EAN of inner box | 4038917409513 |
| EAN case | 4038917409506 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

