



Product description

Order number	77.3898.506
Product description	Disposal container, Multi-Safe medi, 6 l, yellow, round, material: PP, opening: round, height: 195 mm, flat cap, Ø of cap: 230 mm, cap: red, TRBA 250, UN 3291, DIN EN ISO 23907, 21 piece(s)/carton
UN identification	yes

Product characteristics

Cap	flat cap
Shape of opening	round
Removal aid - turning aid	Removal and rotation aid

Size

Max. volume	6 l
Diameter of base	190 mm
Diameter of lid	230 mm
Diameter of opening	88 mm
Height of product	195 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

Purity & certification

Satisfies the requirement	TRBA 250, UN 3291, DIN EN ISO 23907
Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	21
Type of smallest subpackaging	case
Piece(s) / inner box	21
Piece (s) / outer case	21
Piece(s) / pallet	336
Depth of case	395 mm
Width of case	430 mm
Height of case	585 mm
Case volume	0.0994 cbm
Weight of product	0.3572 kg
Weight of case	8 kg
EAN of inner box	4038917543484
EAN case	4038917543491

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

