



### **Product description**

|                     |  |
|---------------------|--|
| Order number        | 77.3899.050  |
| Product description | Disposal container, Multi-Safe eco 50, 50 l, yellow, square, material: PP, opening: square, height: 500 mm, flat cap, lid (LxW): 395 x 295 mm, cap: yellow, UN 3291, 125 piece(s)/carton |
| UN identification   | yes  |

### **Product characteristics**

|                           |               |
|---------------------------|---------------|
| Cap                       | flat cap      |
| Shape of opening          | square        |
| Removal aid - turning aid | not available |

### **Size**

|                     |                 |
|---------------------|-----------------|
| Max. volume         | 50 l            |
| Diameter of base    | 330 mm x 230 mm |
| Diameter of opening | 372 mm x 271 mm |
| Height of product   | 500 mm          |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | yellow             |
| Colour of cap     | yellow             |
| Colour of lid     | yellow             |

### **Purity & certification**

|                           |                         |
|---------------------------|-------------------------|
| Satisfies the requirement | UN 3291                 |
| Product category          | no medical device   IVD |
| Batched                   | no                      |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 125           |
| Type of smallest subpackaging | pallet        |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 125           |
| Piece(s) / pallet             | 125           |
| Case volume                   | 1.72 cbm      |
| Weight of product             | 1.4965 kg     |
| Weight of case                | 187.13 kg     |
| EAN of inner box              | 4038917333528 |
| EAN case                      | 4038917333511 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

