



### **Product description**

Order number	77.3899.060
Product description	Disposal container, Multi-Safe steri, 60 l, yellow, square, material: PP, autoclavable, opening: square, height: 520 mm, flat cap, lid (LxW): 393 x 393 mm, cap: red, UN 3291, 120 piece(s)/carton
UN identification	yes

### **Product characteristics**

Cap	flat cap
Shape of opening	square
Removal aid - turning aid	not available
Autoclavable (max. °C)	134 °C

### **Size**

Max. volume	60 l
Diameter of base	340 mm x 340 mm
Diameter of opening	370 mm x 370 mm
Height of product	520 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	yellow
Colour of cap	red
Colour of lid	yellow

### **Purity & certification**

Satisfies the requirement	UN 3291
Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	120
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	120
Piece(s) / pallet	120
Case volume	2.304 cbm
Weight of product	1.8345 kg
Weight of case	220.2 kg
EAN of inner box	4038917333559
EAN case	4038917333542

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

