



**Product Specification**  
**Post transport packaging, 192 x 146 x 75 mm, for**  
**cold transport container**

**Page 1**



**Product description**

Order number	95.903
Product description	Post transport packaging, application e.g.: external packaging according to P 650/ ADR, (LxWxH): 192 x 146 x 75 mm, for cold transport container, holds 1 to 4 transport containers/bottles up to a diameter of 30 mm, 50 piece(s)/case
UN identification	yes

**Product characteristics**

Label/ Print	with print
Colour of print/label	black
Example of use	external packaging according to P 650/ ADR
Suitable for	cold transport container

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Post transport packaging, 192 x 146 x 75 mm, for**  
**cold transport container**

**Page 2**

**Size**

Width of product	146 mm
Height of product	75 mm
Length of product	192 mm

**Material & colours**

Product material	cardboard
Colour of product	white

**Purity & certification**

Satisfies the requirement	P 650, ADR, IATA
Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	50
Piece(s) / pallet	1650
Depth of case	548 mm
Width of case	513 mm
Height of case	115 mm
Case volume	0.0323 cbm
Weight of product	0.0782 kg
Weight of case	4.54 kg
EAN of inner box	4038917084581
EAN case	4038917079594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-25

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**