



Product description

| | |
|---------------------|----------------------------------|
| Order number | 95.995 |
| Product description | Cold transport container, german |

Product characteristics

| | |
|-----------------------|---------------|
| Label/ Print | with print |
| Colour of print/label | red |
| Application | Cold shipping |
| Absorbent liner | no |

Size

| | |
|-------------------|--------|
| Width of product | 71 mm |
| Height of product | 170 mm |
| Length of product | 140 mm |

Material & colours

| | |
|-------------------|--------------------------------------|
| Product material | Acrylnitrile butadiene styrene (ABS) |
| Colour of product | white |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | P 650, ADR |
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 6 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 6 |
| Piece(s) / pallet | 648 |
| Depth of case | 308 mm |
| Width of case | 223 mm |
| Height of case | 195 mm |
| Case volume | 0.0134 cbm |
| Weight of product | 0.3448 kg |
| Weight of case | 2.3 kg |
| EAN of inner box | 4038917084901 |
| EAN case | 4038917003667 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

