



Packaging

Minimum order qty.	2000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	80000
Depth of case	428 mm
Width of case	388 mm
Height of case	225 mm
Case volume	0.0374 cbm
Weight of product	0.0024 kg
Weight of case	5.47 kg
EAN of inner box	4038917082723
EAN case	4038917001762

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

