



Product description

Order number	84.1412.033
Product description	Adhesive label, (LxW): 145 x 98.5 mm, white, "EXEMPT MEDICAL SAMPLE", print: black, removable, 1,000 piece(s)/roll

Product characteristics

Colour of print/label	black
-----------------------	-------

Size

Width of product	98.5 mm
Length of product	145 mm

Material & colours

Colour of product	white
-------------------	-------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest roll
subpackaging
Piece(s) / inner box 1000
Piece (s) / outer case 1000
Piece(s) / pallet 144000
Depth of case 323 mm
Width of case 213 mm
Height of case 145 mm
Case volume 0.01 cbm
Weight of product 0.0022 kg
Weight of case 2.5 kg
EAN of inner box 4038917223119
EAN case 4038917223102

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

