



Product description

Order number	86.9925.424
Product description	Specimen bag, (LxW): 260 x 160 mm, material: LD-PE, transparent, film thickness: 40 µm, print Biohazard, 100 piece(s)/bag

Product characteristics

Application	Transportation of diagnostic specimens
-------------	--

Size

Width of product	160 mm
Length of product	260 mm

Material & colours

Product material	Low Density Polyethylene (LD-PE)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Specimen bag, (LxW): 260 x 160 mm, LD-PE,
transparent, Biohazard

Purity & certification

Latex-free	yes
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	96000
Depth of case	428 mm
Width of case	358 mm
Height of case	234 mm
Case volume	0.0357 cbm
Weight of case	13.3 kg
EAN of inner box	4038917155557
EAN case	4038917040129

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

