



### **Product description**

Order number	95.978
Product description	Mailing envelope, material: PE, (LxW): 230 x 300 mm, with 50 mm adhesive tape, film thickness: 0.045 mm, 100 piece(s)/case

### **Product characteristics**

Design	with 50 mm adhesive tape
Absorbent liner	no

### **Size**

Width of product	300 mm
Length of product	230 mm

### **Material & colours**

Product material	Polyethylene (PE)
------------------	-------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	56000
Depth of case	290 mm
Width of case	145 mm
Height of case	400 mm
Case volume	0.0168 cbm
Weight of product	0.0057 kg
Weight of case	6.12 kg
EAN of inner box	4038917224178
EAN case	4038917224161

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

