




Product Specification
Superabsorber pack, suitable for transport case
T 15 and B 17, (LxW): 400 x 300 mm

Page 1



Product description

Order number	95.1712
Product description	Superabsorber pack, suitable for transport case T 15 and B 17, (LxW): 400 x 300 mm, fluid intake: 1.2 l, 100 piece(s)/case

Product characteristics

Application	Mailing of samples, Multi-Safe disposel boxes
Suitable for	transport case T 15 and B 17

Size

Width of product	300 mm
Length of product	400 mm

Material & colours

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Superabsorber pack, suitable for transport case
T 15 and B 17, (LxW): 400 x 300 mm

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	100
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	100
Piece(s) / pallet	72000
Depth of case	420 mm
Width of case	320 mm
Height of case	140 mm
Case volume	0.0188 cbm
Weight of product	0.0283 kg
Weight of case	3.5 kg
EAN of inner box	4038917415774
EAN case	4038917065894

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

