



### **Product description**

Order number	95.1714
Product description	Closure clip, (LxW): 220 x 12 mm, material: PA, for cross-bottom bag 95.1713, 25 piece(s)/case

### **Product characteristics**

Suitable for	cross-bottom bag 95.1713
--------------	--------------------------

### **Size**

Width of product	12 mm
Length of product	220 mm

### **Material & colours**

Product material	Polyamide (PA)
------------------	----------------

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Packaging**

Minimum order qty.	25
Type of smallest subpackaging	case
Piece(s) / inner box	25
Piece (s) / outer case	25
Piece(s) / pallet	8100
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0512 kg
Weight of case	1.4 kg
EAN of inner box	4038917087520
EAN case	4038917065948

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

