



**Product Specification**  
**Transport case B 17, with lateral rim for block rack**

**Page 1**



**Product description**

Order number	95.1717
Product description	Transport case B 17, with lateral rim for block rack, (LxWxH): 395 x 395 x 435 mm

**Product characteristics**

Design	with lateral rim for block rack
--------	---------------------------------

**Size**

Width of product	395 mm
Height of product	435 mm
Length of product	395 mm

**Material & colours**

Product material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Transport case B 17, with lateral rim for block rack**

**Purity & certification**

Satisfies the requirement	P 650, ADR
Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	18
Depth of case	403 mm
Width of case	403 mm
Height of case	460 mm
Case volume	0.0747 cbm
Weight of product	3.561 kg
Weight of case	4.5 kg
EAN of inner box	4038917087544
EAN case	4038917065610

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

