



Product description

Order number	95.1725
Product description	Battery bag, without batteries, (LxW): 205 x 115 mm, material: PL, for standard cooling batteries, 12 piece(s)/case

Product characteristics

Design	without batteries
Suitable for	standard cooling batteries

Size

Width of product	115 mm
Length of product	205 mm

Material & colours

Product material	polyester (PL)
------------------	----------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 12
Type of smallest case
subpackaging
Piece(s) / inner box 12
Piece (s) / outer case 12
Piece(s) / pallet 1980
Depth of case 338 mm
Width of case 228 mm
Height of case 120 mm
Case volume 0.0092 cbm
Weight of product 0.4419 kg
Weight of case 5.5 kg
EAN of inner box 4038917087568
EAN case 4038917068376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

