



**Product Specification**  
**Screw cap, natural, suitable for mailing container**  
**85 x 30 mm**

**Page 1**



**Product description**

Order number	65.676
Product description	Screw cap, natural, suitable for mailing container 85 x 30 mm, 50 piece(s)/bag

**Product characteristics**

Application	Mailing of samples
Suitable for	mailing container 85 x 30 mm

**Size**

Length of product	15 mm
-------------------	-------

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	natural

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap, natural, suitable for mailing container**  
**85 x 30 mm**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      500  
Type of smallest                              bag  
subpackaging  
Piece(s) / inner box                      50  
Piece (s) / outer case                      500  
Piece(s) / pallet                              39000  
Depth of case                              418 mm  
Width of case                                278 mm  
Height of case                                130 mm  
Case volume                                0.0151 cbm  
Weight of product                        0.0028 kg  
Weight of case                                1.73 kg  
EAN of inner box                            4038917083300  
EAN case                                      4038917002240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

