




Product Specification
Screw cap, brown, suitable for Mailing container
126 x 30 mm

Page 1



Product description

| | |
|---------------------|---|
| Order number | 65.679.001 |
| Product description | Screw cap, brown, suitable for Mailing container 126 x 30 mm, 50 piece(s)/bag |

Product characteristics

| | |
|--------------|-------------------------------|
| Application | Mailing of samples |
| Suitable for | Mailing container 126 x 30 mm |

Size

| | |
|-------------------|-------|
| Length of product | 20 mm |
|-------------------|-------|

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | brown |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Screw cap, brown, suitable for Mailing container
126 x 30 mm

Page 2

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 250 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 250 |
| Piece(s) / pallet | 41250 |
| Depth of case | 338 mm |
| Width of case | 228 mm |
| Height of case | 120 mm |
| Case volume | 0.0092 cbm |
| Weight of product | 0.0034 kg |
| Weight of case | 1.09 kg |
| EAN of inner box | 4038917115810 |
| EAN case | 4038917014915 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT