




Product Specification
**Absorbent liner, suitable for Mailing container 114
x 44 mm, (LxW): 135 x 95 mm**

Page 1



Product description

Order number	78.669
Product description	Absorbent liner, suitable for Mailing container 114 x 44 mm, (LxW): 135 x 95 mm, 100 piece(s)/case

Product characteristics

Application	Mailing of samples, Multi-Safe disposel boxes
Suitable for	Mailing container 114 x 44 mm

Size

Width of product	95 mm
Length of product	135 mm

Material & colours

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Absorbent liner, suitable for Mailing container 114
x 44 mm, (LxW): 135 x 95 mm

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1000
Type of smallest case
subpackaging
Piece(s) / inner box 100
Piece (s) / outer case 1000
Piece(s) / pallet 56000
Depth of case 385 mm
Width of case 385 mm
Height of case 160 mm
Case volume 0.0237 cbm
Weight of product 0.0026 kg
Weight of case 3.04 kg
EAN of inner box 4038917083249
EAN case 4038917002196

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

