



**Product Specification**  
**Absorbent liner, suitable for Mailing container 126  
x 30 mm, (LxW): 75 x 90 mm**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 78.670   |
| Product description | Absorbent liner, suitable for Mailing container 126 x 30 mm,<br>(LxW): 75 x 90 mm, 1,000 piece(s)/case |

**Product characteristics**

|              |   |
|--------------|---|
| Application  | Mailing of samples, Multi-Safe disposel boxes |
| Suitable for | Mailing container 126 x 30 mm                 |

**Size**

|                   |       |
|-------------------|-------|
| Width of product  | 90 mm |
| Length of product | 75 mm |

**Material & colours**

|                  |           |
|------------------|-----------|
| Product material | cellulose |
|------------------|-----------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Absorbent liner, suitable for Mailing container 126**  
**x 30 mm, (LxW): 75 x 90 mm**

**Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 4000          |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 1000          |
| Piece (s) / outer case        | 4000          |
| Piece(s) / pallet             | 72000         |
| Depth of case                 | 515 mm        |
| Width of case                 | 360 mm        |
| Height of case                | 360 mm        |
| Case volume                   | 0.0667 cbm    |
| Weight of product             | 0.001 kg      |
| Weight of case                | 4.36 kg       |
| EAN of inner box              | 4038917377492 |
| EAN case                      | 4038917159937 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

