




**Product Specification**  
**Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm**

**Page 1**



**Product description**

Order number	78.671
Product description	Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm, 100 piece(s)/case

**Product characteristics**

Application	Mailing of samples, Multi-Safe disposel boxes
Suitable for	cold transport container

**Size**

Width of product	55 mm
Length of product	80 mm

**Material & colours**

Product material	cellulose
------------------	-----------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Absorbent liner, suitable for cold transport container, (LxW): 80 x 55 mm**

**Page 2**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

**Packaging**

Minimum order qty.                      2500  
Type of smallest  
subpackaging                              case  
Piece(s) / inner box                      100  
Piece (s) / outer case                      2500  
Piece(s) / pallet                           270000  
Depth of case                              383 mm  
Width of case                               193 mm  
Height of case                              190 mm  
Case volume                                0.014 cbm  
Weight of product                        0.0009 kg  
Weight of case                              2.57 kg  
EAN of inner box                          4038917083263  
EAN case                                    4038917079860

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**