



Product Specification

**Mailing container, transparent, construction:
round, length: 114 mm, Ø opening: 44 mm,
without cap**

Page 1



Product description

| | |
|---------------------|---|
| Order number | 78.895 |
| Product description | Mailing container, round, transparent, construction: round, length: 114 mm, Ø opening: 44 mm, without cap, cap: natural, 50 piece(s)/case |

Product characteristics

| | |
|------------------|--|
| Application | Mailing of samples |
| Example of use | Secondary packaging according to P 650 |
| Absorbent liner | no |
| Shape of opening | round |
| Liquid-tight | yes |

Size

| | |
|---------------------|--------|
| Diameter of opening | 44 mm |
| Length of product | 114 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**Mailing container, transparent, construction:
round, length: 114 mm, Ø opening: 44 mm,
without cap**

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Colour of cap | natural |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | P 650, ADR |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 250 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 250 |
| Piece(s) / pallet | 5000 |
| Depth of case | 598 mm |
| Width of case | 398 mm |
| Height of case | 350 mm |
| Case volume | 0.0833 cbm |
| Weight of product | 0.017 kg |
| Weight of case | 5.25 kg |
| EAN of inner box | 4038917084512 |
| EAN case | 4038917003322 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

