



Product Specification
**Urine container, 3 l, with integrated transfer unit,
brown, with light protection, graduated**

Page 1



Product description

Order number	77.575.402
Product description	Urine container, Urine collection, max. working volume: 3 l, Ø opening: 83.5 mm, brown, with light protection, graduated, material: PE, with plastic label, French, cap: yellow, with integrated transfer unit, 30 piece(s)/case

Product characteristics

Label/ Print	with plastic label
Graduation	yes
Closure type	screw cap
Light protection	yes
Design	with integrated transfer unit
Application	Urine collection
Example of use	closed urine transfer into a V-Monovette® Urine
Volume of work	700 ml
Viewing window	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine container, 3 l, with integrated transfer unit,
brown, with light protection, graduated**

Page 2

Size

Sample volume	3.000.00 ml
Max. volume	3 l
Diameter of opening	83,5 mm
Width of product	130 mm
Length of product	125 mm
Length including cap	246 mm
Length excluding cap	243 mm

Material & colours

Product material	Polyethylene (PE)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	yellow

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	30
Type of smallest subpackaging	case
Piece(s) / inner box	30
Piece (s) / outer case	30
Piece(s) / pallet	270
Depth of box	643 mm
Width of box	403 mm
Height of box	515 mm
Depth of case	643 mm
Width of case	403 mm
Height of case	515 mm
Case volume	0.1335 cbm
Weight of product	0.1433 kg
Weight of case	5.1 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine container, 3 l, with integrated transfer unit,
brown, with light protection, graduated**

Page 3

EAN of inner box	4038917344791
EAN case	4038917344784

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT