



**Product Specification**  
**Urine container, 2 l, brown, with light protection,  
graduated**

**Page 1**



**Product description**

Order number	77.580
Product description	Urine container, Urine collection, max. working volume: 2 l, Ø opening: 70 mm, brown, with light protection, graduated, material: PE, with paper label, german/english, cap: yellow, Label loosely enclosed, 30 piece(s)/case

**Product characteristics**

Label/ Print	with paper label
Colour of print/label	yellow/black
Graduation	yes
Closure type	screw cap
Light protection	yes
Application	Urine collection
Viewing window	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine container, 2 l, brown, with light protection,  
graduated**

**Page 2**

**Size**

Sample volume	2.000.00 ml
Max. volume	2 l
Diameter of opening	70 mm
Width of product	106 mm
Height of product	240 mm
Length of product	106 mm
Length including cap	238 mm
Length excluding cap	235 mm

**Material & colours**

Product material	Polyethylene (PE)
Colour of product	brown
Closure material	Polypropylene (PP)
Colour of cap	yellow

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Safety data sheet	no
Batched	yes

**Packaging**

Minimum order qty.	30
Type of smallest subpackaging	case
Piece(s) / inner box	30
Piece (s) / outer case	30
Piece(s) / pallet	360
Depth of case	563 mm
Width of case	343 mm
Height of case	500 mm
Case volume	0.0966 cbm
Weight of product	0.0956 kg
Weight of case	3.5 kg
EAN of inner box	4038917082655
EAN case	4038917001700

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

