



## Product Specification

Medicine cup, max. volume: 30 ml, (LxØ): 40 x 37 mm, graduated, PP, yellow

Page 1



### Product description

Order number	75.1337.504
Product description	Medicine cup, max. volume: 30 ml, (LxØ): 40 x 37 mm, Ø opening: 37 mm, yellow, graduated, material: PP, 75 piece(s)/bag

### Product characteristics

Graduation	yes
Volume of work	30 ml

### Size

Sample volume	30 ml
Max. volume	30 ml
Diameter	37 mm
Diameter of opening	37 mm
Height of product	40 mm
Length of product	40 mm
Length excluding cap	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Medicine cup, max. volume: 30 ml, (LxØ): 40 x 37 mm, graduated, PP, yellow**

**Page 2**

### Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow

### Purity & certification

Product category	Medical device (MD)
CE certified	CE supplier
Batched	yes

### Packaging

Minimum order qty.	3750
Type of smallest subpackaging	bag
Piece(s) / inner box	3750
Piece (s) / outer case	3750
Piece(s) / pallet	180000
Depth of box	190 mm
Width of box	390 mm
Height of box	390 mm
Depth of case	400 mm
Width of case	400 mm
Height of case	205 mm
Case volume	0.0328 cbm
Weight of product	0.0015 kg
Weight of case	6 kg
EAN of inner box	4038917127745
EAN case	4042836002052

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**