



Product description

Order number	75.1354.001
Product description	Container with screw cap, collection and storage of urine, max. working volume: 100 ml, Ø opening: 58 mm, transparent, graduated, material: PP, screw cap, cap: green, cap enclosed, 25 pieces each stacked inside each other, 500 piece(s)/bag

Product characteristics

Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Volume of work	100 ml

Size

Sample volume	100 ml
Max. volume	100 ml
Diameter of opening	58 mm
Length including cap	73 mm
Length excluding cap	71 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE supplier
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	500
Piece(s) / pallet	12000
Depth of case	600 mm
Width of case	330 mm
Height of case	400 mm
Case volume	0.0778 cbm
Weight of product	0.0115 kg
Weight of case	7.64 kg
EAN of inner box	4038917247801
EAN case	8052109529632

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

