



Product description

Order number	75.1355
Product description	Urine container, collection and storage of urine, max. working volume: 125 ml, (ØxH): 66 x 67 mm, Ø opening: 66 mm, highly transparent, graduated, material: PP, 50 piece(s)/bag

Product characteristics

Graduation	yes
Application	collection and storage of urine
Volume of work	125 ml

Size

Sample volume	125 ml
Max. volume	125 ml
Diameter	66 mm
Diameter of opening	66 mm
Height of product	67 mm
Length excluding cap	67 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Urine container, 125 ml, (ØxH): 66 x 67 mm, PP,
highly transparent**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	highly transparent

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE supplier
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	50
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of box	500 mm
Width of box	65 mm
Height of box	65 mm
Depth of case	590 mm
Width of case	400 mm
Height of case	280 mm
Case volume	0.0661 cbm
Weight of product	0.0057 kg
Weight of case	6.2 kg
EAN of inner box	4038917409568
EAN case	4042836010088

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT