




**Product Specification**  
**Urine container, 90 ml, (ØxH): 60 x 65 mm, PP,**  
**with paper label, transparent**



**Product description**

Order number	75.1435.004
Product description	Urine container, max. working volume: 90 ml, (ØxH): 60 x 65 mm, Ø opening: 65 mm, transparent, graduated, material: PP, screw cap, with paper label, cap assembled, sterile, 400 piece(s)/case

**Product characteristics**

Label/ Print	with paper label
Graduation	yes
Closure type	screw cap
Volume of work	90 ml

**Size**

Sample volume	90 ml
Max. volume	90 ml
Diameter	60 mm
Diameter of opening	65 mm
Height of product	65 mm
Length excluding cap	60 µm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Urine container, 90 ml, (ØxH): 60 x 65 mm, PP,  
with paper label, transparent**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	400
Piece (s) / outer case	400
Piece(s) / pallet	7200
Depth of case	610 mm
Width of case	286 mm
Height of case	610 mm
Case volume	0.1062 cbm
Weight of product	0.0201 kg
Weight of case	8.23 kg
EAN of inner box	4038917385220
EAN case	4038917385213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**