



Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated transfer unit, transparent

Page 1



Product description

| | |
|---------------------|--|
| Order number | 75.562.400 |
| Product description | Container with screw cap, collection and storage of urine, max. working volume: 100 ml, (ØxH): 57 x 76 mm, Ø opening: 62 mm, transparent, graduated, material: PP, screw cap, cap: yellow, with safety label, with integrated transfer unit, cap assembled, sterile, 5 piece(s)/tube bag |

Product characteristics

| | |
|-----------------------|---|
| Label/ Print | with safety label |
| Colour of print/label | white/red/black |
| Graduation | yes |
| Closure type | screw cap |
| Design | with integrated transfer unit |
| Application | collection and storage of urine |
| Example of use | closed urine transfer into a V-Monovette® Urine |
| Volume of work | 20 ml - 100 ml |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Container with screw cap, 100 ml, (ØxH): 57 x 76 mm, PP, with safety label, with integrated transfer unit, transparent

Page 2

Size

| | |
|----------------------|--------|
| Sample volume | 100 ml |
| Max. volume | 100 ml |
| Diameter | 57 mm |
| Diameter of opening | 62 mm |
| Height of product | 76 mm |
| Length excluding cap | 72 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | yellow |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 200 |
| Type of smallest subpackaging | tube bag |
| Piece(s) / inner box | 200 |
| Piece (s) / outer case | 200 |
| Piece(s) / pallet | 5600 |
| Depth of case | 588 mm |
| Width of case | 398 mm |
| Height of case | 255 mm |
| Case volume | 0.0597 cbm |
| Weight of product | 0.0173 kg |
| Weight of case | 4.31 kg |
| EAN of inner box | 4038917236805 |
| EAN case | 4038917236799 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

