



Product Specification
Multi-purpose container, max. volume: 26 ml,
(LxØ): 45 x 70 mm, PP, transparent

Page 1



Product description

Order number	75.568
Product description	Multi-purpose container, max. volume: 26 ml, (LxØ): 45 x 70 mm, Ø opening: 70 mm, transparent, material: PP

Product characteristics

Graduation	no
Application	collection and storage of urine
Volume of work	26 ml

Size

Sample volume	26 ml
Max. volume	26 ml
Diameter	70 mm
Diameter of opening	70 mm
Height of product	45 mm
Length of product	45 mm
Length excluding cap	45 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Multi-purpose container, max. volume: 26 ml,
(LxØ): 45 x 70 mm, PP, transparent

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category not regulated
CE certified no CE mark
Batched no

Packaging

Minimum order qty. 1000
Piece(s) / inner box 500
Piece (s) / outer case 1000
Piece(s) / pallet 30000
Depth of case 398 mm
Width of case 398 mm
Height of case 340 mm
Case volume 0.0539 cbm
Weight of product 0.0036 kg
Weight of case 4.25 kg
EAN of inner box 4038917082549
EAN case 4038917001625

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT