




## Product Specification

**Multi-purpose container, max. volume: 70 ml,  
(LxØ): 55 x 44 mm, graduated, PP, transparent,  
with paper label**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 75.9922.721  |
| Product description | Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, with paper label, sterile, 1 piece(s)/blister |

### Product characteristics

|                |                                 |
|----------------|---------------------------------|
| Label/ Print   | with paper label                |
| Graduation     | yes                             |
| Closure type   | screw cap                       |
| Application    | collection and storage of urine |
| Base shape     | flat base                       |
| Volume of work | 50 ml                           |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Multi-purpose container, max. volume: 70 ml,  
(LxØ): 55 x 44 mm, graduated, PP, transparent,  
with paper label**

**Page 2**

### Size

|                      |       |
|----------------------|-------|
| Max. volume          | 70 ml |
| Diameter             | 44 mm |
| Diameter of opening  | 44 mm |
| Height of product    | 55 mm |
| Length of product    | 55 mm |
| Length excluding cap | 55 mm |

### Material & colours

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polypropylene (PP)                |
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | yellow                            |

### Purity & certification

|                  |                                      |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD)            |
| CE certified     | CE - manufacturer's self-declaration |
| Purity standard  | sterile                              |
| Sterilisation    | Gamma irradiation                    |
| Batched          | yes                                  |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 200           |
| Type of smallest subpackaging | blister       |
| Piece(s) / inner box          | 200           |
| Piece (s) / outer case        | 200           |
| Piece(s) / pallet             | 3200          |
| Depth of case                 | 548 mm        |
| Width of case                 | 375 mm        |
| Height of case                | 405 mm        |
| Case volume                   | 0.0832 cbm    |
| Weight of product             | 0.0112 kg     |
| Weight of case                | 3.5 kg        |
| EAN of inner box              | 4038917152525 |
| EAN case                      | 4038917347181 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

