



## Product Specification

**Multi-purpose container, max. volume: 70 ml,  
(LxØ): 55 x 44 mm, graduated, PP, transparent**

**Page 1**



### Product description

Order number	75.9922.744
Product description	Multi-purpose container, max. volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, 250 piece(s)/bag

### Product characteristics

Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Base shape	flat base
Volume of work	50 ml

### Size

Max. volume	70 ml
Diameter	44 mm
Diameter of opening	44 mm
Height of product	55 mm
Length of product	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**Multi-purpose container, max. volume: 70 ml,  
(LxØ): 55 x 44 mm, graduated, PP, transparent**

**Page 2**

Length excluding cap 55 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	8000
Depth of case	555 mm
Width of case	380 mm
Height of case	405 mm
Case volume	0.0854 cbm
Weight of product	0.0123 kg
Weight of case	7 kg
EAN of inner box	4038917152754
EAN case	4038917027793

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**