



**Product Specification**  
**Multi-purpose container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent**

**Page 1**



**Product description**

Order number	75.9922.744
Product description	Multi-purpose container, max. working volume: 70 ml, (LxØ): 55 x 44 mm, Ø opening: 44 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, 250 piece(s)/bag

**Product characteristics**

Graduation	yes
Closure type	screw cap
Application	collection and storage of urine
Base shape	flat base
Volume of work	70 ml

**Size**

Sample volume	70 ml
Max. volume	70 ml
Diameter	44 mm
Diameter of opening	44 mm
Height of product	55 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Multi-purpose container, 70 ml, (LxØ): 55 x 44 mm, graduated, PP, transparent**

**Page 2**

Length of product 55 mm  
Length excluding cap 55 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap yellow

**Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest subpackaging bag  
Piece(s) / inner box 250  
Piece (s) / outer case 500  
Piece(s) / pallet 8000  
Depth of case 548 mm  
Width of case 408 mm  
Height of case 373 mm  
Case volume 0.0834 cbm  
Weight of product 0.0123 kg  
Weight of case 7 kg  
EAN of inner box 4038917152754  
EAN case 4038917027793

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**