



Product description

| | |
|---------------------|--|
| Order number | 76.1340.500 |
| Product description | Lid, for 30-ml medicine cup, material: PE, white, 750 piece(s)/bag |

Product characteristics

| | |
|----------------|---|
| Example of use | Oral administration of medications in liquid or tablet form |
| Suitable for | 30-ml medicine cup |

Material & colours

| | |
|-------------------|-------------------|
| Product material | Polyethylene (PE) |
| Colour of product | white |
| Colour of cap | white |

Purity & certification

| | |
|------------------|----------------|
| Product category | Medical device |
| CE certified | CE supplier |
| Batched | yes |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 3750 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 3750 |
| Piece (s) / outer case | 3750 |
| Piece(s) / pallet | 180000 |
| Depth of case | 190 mm |
| Width of case | 390 mm |
| Height of case | 390 mm |
| Case volume | 0.0289 cbm |
| Weight of product | 0.0007 kg |
| Weight of case | 2.97 kg |
| EAN of inner box | 4038917127790 |
| EAN case | 4042836007149 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

