



Product Specification
Snap-on lid, with spout, PE, red, 1,000
piece(s)/bag

Page 1



Product description

Order number	76.1477
Product description	Snap-on lid, with spout, suitable for 125-ml container 75.1355, material: PE, red, 1,000 piece(s)/bag

Product characteristics

Closure type	snap-on lid
Design	with spout
Suitable for	125-ml container 75.1355

Material & colours

Product material	Polyethylene (PE)
Colour of product	red

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE supplier
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Snap-on lid, with spout, PE, red, 1,000
piece(s)/bag

Page 2

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	1000
Piece(s) / pallet	24000
Depth of case	390 mm
Width of case	265 mm
Height of case	590 mm
Case volume	0.061 cbm
Weight of product	0.0037 kg
Weight of case	4.47 kg
EAN of inner box	4038917086752
EAN case	4042836001574

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT