



### **Product description**

Order number	57.462
Product description	Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Closure type	push cap
Base shape	conical base
Volume of work	12 ml
Centrifugation max (RCF)	2500 x g

### **Size**

Diameter	17 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material                      Polystyrene (PS)  
Colour of product                      transparent

### **Purity & certification**

Satisfies the requirement              sealed at 0.6 bar/10 min  
Product category                      In-vitro diagnostic (IVD)  
CE certified                              CE - manufacturer's self-declaration  
Batched                                      yes

### **Packaging**

Minimum order qty.                      1000  
Type of smallest                          bag  
subpackaging  
Piece(s) / inner box                      500  
Piece (s) / outer case                    1000  
Piece(s) / pallet                          36000  
Depth of case                              398 mm  
Width of case                               398 mm  
Height of case                              310 mm  
Case volume                               0.0491 cbm  
Weight of product                        0.0063 kg  
Weight of case                              6.98 kg  
EAN of inner box                          4038917081702  
EAN case                                    4038917000819

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

